# BIOpick

*Your Potato Sample Device*

## USER GUIDE

## Table of contents

1	Foreword	3
2	Content check	3
3	Safety warning and guidelines	3
4	Decontamination of BIOpick and collection tubes	3
5	Introduction	4
Table 1.	Specifications	4
6	Operation guide	4
6.1	BIOpick: manual use	4
6.2	BIOpick: semi-automated use	4
7	Frequently Asked Questions	5
8	Video tutorials	5
9	Ordering information	5

# 1. Foreword

Thank you for purchasing the BIOpick. The BIOpick is a device for collecting heel core samples from potato tubers for downstream analysis of potato pathogens.

It makes coring more ergonomic allowing collection and transfer of multiple samples into a sterile 50 ml collection tube (Fig. 1) in a single step. The collection of core samples can be done manually or semi-automatically by connecting the BIOpick to a laboratory mixer through the adapter for BIOpick (Fig. 2).

# 2. Content check

Please check the content (see chapter 9) when you first open the BIOpick packaging. If items are damaged or missing, please contact us immediately.


# 3. Safety warnings and guidelines

The operation and maintenance of this instrument should comply with this user guide. Failure to comply with these guidelines may affect the durability of the device. Do not use the BIOpick and the adapter for anything other than its intended use.

# 4. Decontamination of BIOpick and collection tubes

To avoid cross-contaminations, the BIOpick must be washed thoroughly in a cycle of 0.5% sodium hypochlorite (household bleach) for about 30 minutes, 70% ethanol, and distilled water\* or according to your local guidelines. The BIOpick can also be autoclaved at 121 °C for 30-60 min. Sample leftovers inside of the screw cap should be removed using the provided cleaning brush. When not used, the BIOpick and accessories should be stored in a dry place. Some polypropylene 50 ml screw-cap tubes can be autoclaved at 121 °C (not the caps) without any appreciable loss of mechanical properties and can therefore be re-used for sample collection\*\*.

\*Our tests showed that after treatment with household bleach, ethanol and dest. water in combination with a cleaning brush, Potato Virus Y (PVY) was not detectable by qPCR when healthy heel samples were collected with a decontaminated BIOpick previously used to collect PVY-infected samples (see graph below).

\*\*Our tests showed that it was possible to autoclave (121°C for 30 min) both BIOpick and screw cap tubes of the brand Falcon (VWR, Art. Nr. 734-0448) multiple times without impairing the functionality of either of the BIOpick nor the tubes.

## 5. Introduction

The BIOpick (Table 1) has been designed specifically to collect potato heel cores samples in the context of seed potatoes certification according to EU standards (commission implementing regulation 2022/1194 of 11 July 2022). To detect quarantine bacterial pests such as potato ring rot (Cms; *Clavibacter michiganensis* subsp. *sepedonicus*) and potato brown rot (Rs; *Ralstonia solanacearum*), EU legislation defines a standard sample size of 200 heel-end cores. The BIOpick can also be used to collect 6 x 50 heel-end cores for the detection of *Dickeya* spp. and *Pectobacteria* spp. causing blackleg or for virus detection. Connecting the BIOpick to a laboratory mixer significantly reduces the effort and the time of sampling (watch a video in section 8).

Trademark	BIOpick
Model	Potato Sampler
Size of core	4 mm (diameter) x 5-7 mm (length)
Weight of core	~ 100 mg per core
Dimensions	BIOpick: 25.5 mm (height) x 36 mm (diameter) Adapter: 130 mm (length) x 34 mm (diameter)
Weight	BIOpick: 40 g Adapter: 323.8 g
Material	Stainless steel

Table 1. Specifications

## 6. Operation guide

The BIOpick can be used manually (see 6.1) or semi-automatically (see 6.2) in combination with a laboratory mixer with adjustable rotation speed (not provided). We recommend both laboratory mixer (RSO-20A) and tripod (RSO-E07) from Phoenix Instrument GmbH but other models can be used as well. The following paragraphs (6.1. and 6.2) describe the recommended steps to collect heel cores for Cms/Rs testing, however they can be modified depending on the pathogen to be tested.

### 6.1 BIOpick: manual use

- Take a screw cap tube (50 ml) and replace the plastic cap with the BIOpick (Fig. 3).
- Remove the epidermis around the heel from a potato tuber using a peeler (Fig. 4). You can peel the tubers ahead of time to make the process smoother.
- Insert the BIOpick into the heel tissue to collect a heel core and repeat this step with the next potato tuber (Fig. 5).
- Collect up to 200 heel cores in one 50 ml screw cap tube. When switching to different samples, it is important to disinfect all parts to avoid cross-contaminations (see chapter 4).

Watch a video-tutorial on how to collect samples manually in chapter 8.

### 6.2 BIOpick: semi-automated use

- Take a screw cap tube (50 ml) and replace the plastic cap with the BIOpick (Fig. 3).
- Insert the tube with the BIOpick into the provided adapter and secure it by turning the screw as shown in Fig. 6.
- Attach the adapter to a laboratory mixer (not provided), as shown in Fig. 2.

- Set the speed according to your needs (e.g. 350 rpm). Rotation has to be counterclockwise.
- Remove the epidermis of the tuber with the help of a peeler (you can peel the tubers ahead of time to make the process smoother) (Fig. 4).
- Push the heel end towards the tip of the rotating device. The BIOpick extracts the heel tissue from the tuber (Fig. 7).
- Proceed with the collection of up to 200 heels in one 50 ml screw cap tube. When switching to different samples, it is important to disinfect all parts to avoid cross-contaminations (see chapter 4).

Watch a video-tutorial on how to collect samples semi-automatically in chapter 8.

## 7. Frequently Asked Questions

Problem / Question	Possible cause	Solution
The device does not extract the heel tissue from the tuber when used semi-automatically	The rotation speed of the laboratory mixer may be too high	Slow down the speed of the laboratory mixer and push the tuber more slowly against the tip
The last collected heel core is stuck in the BIOpick		Use a sterile tip or a toothpick to push the heel core into the screw cap tube

## 8. Video tutorials

BIOpick: Manual use

BIOpick: Semi-automatic use

## 9. Ordering information

ITEM	AMOUNT
• BIOpick set (BIOpick, Adapter, Cleaning brushes)(Art. No. BPS01)	1 pc
• BIOpick (Art. No. BP003)	3 pcs
• Adapter for BIOpick (Art. No. BPA01)	1 pc
• Cleaning brushes (Art. No. BPB03)	3 pcs

Fig. 1


Fig. 2

Fig. 3

Fig. 4


Fig. 5

Fig. 6

Fig. 7

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